



# Panther® Link\*

Virtually connects multiple Panther® systems to share information and monitor from a centralised dashboard.

**Value to your lab.** Includes Panther® system savings plus the following incremental savings:

  
LABOUR SAVINGS

Monitor inventories, maintenance and messages from a **centralised dashboard** reducing visits to each system

**Review** reports from a centralised location, maximising tech time

**Share master lot, EQC and user account** information across multiple instruments, to minimise manual entry

**Track and trend QC data for all instruments** to support internal quality procedures and metrics

  
TIME SAVINGS

**Share reflex test orders** (i.e., HPV 16 18/45) across all instruments, providing flexible options for sample loading

Utilise **remote diagnostic** capabilities to help ensure faster troubleshooting and instrument uptime

**Share assay reagent kits** between instruments for better utilisation and to reduce preparation and storage

  
COST SAVINGS

Prevent **inadvertent reloading** of samples with valid test results, avoiding stressful operator errors

Leverage existing Panther investment with a **simple in-lab upgrade**



Adds communication between Panther instruments to share inventories, reagents and results.

Specifications

Dashboard server

Processor	Intel Core i5-6500TE (Skylake) 2.3 GHz processor
Memory	16 GB SO-DIMM DDR4 2133
Primary storage	1 TB 2.5" SSD
Network	Dual 1 GB ethernet
Security	AppArmor
Operating system	Ubuntu 18.04 LTS
Dimensions	19.6 cm W x 5.7 cm D x 26.9 cm H

Dashboard client PC


Processor	Intel Braswell Celeron N3160
Memory	4 GB SO-DIMM DDR3L 1600
Primary storage	32 GB mSATA SSD
Network	1 GB ethernet
Operating system	Porteus Kiosk 4.8.0
Dimensions	11.7 cm W x 3.7 cm D x 11.2 cm H
Mounting brackets	Wall mount

Dashboard command display

Browsers supported	Chrome or mobile Safari (iPad only)
Max. concurrent web clients	50
Resolution	4K (3840 x 2160)
Model/dimensions	QM49R Samsung 49" or QM55R Samsung 55" or QM65R Samsung 65" or Dell P2415Q 23.8" 4K UHD and options to access through an iPad (iOS 12 or later/2018 iPad or newer)

Firewall

Model	Cisco ASA 5506-X
Dimensions	20 cm W x 4.4 cm D x 23.4 cm H
Memory	4 GB
Primary storage	50 GB mSata SSD



Access instrument status and messages from an iPad.



Monitor instrument information and results from an office computer.